
DESS

> DENTAL SMART SOLUTIONS

IFU

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Prosthetic Fixtures for dental implants

INSTRUCTIONS FOR USE (INTERNATIONAL)

ENGLISH

CAUTION – PLEASE READ CAREFULLY

DESS® products should only be used by dental specialists with experience in maxillary implantology and other specialties, such as dental diagnosis, planning, dental surgery, or prosthetic techniques. If in doubt regarding the product's use, please contact the manufacturer. DESS® products designed for single use only must never be re-used. If reused, there is a risk that product damage and deterioration of its characteristics could lead to prosthetic solution failure and / or other deterioration of the patient's health such as patient tissue infection. All DESS® components must be dry fitted before use to check that the fit is correct. The clinician is responsible for the correct application of restorative products as planning and procedures are under his/her control. This is the reason why only dental specialists with the appropriate experience and training should work with these products. In case of doubt please contact the manufacturer or your local distributor.

U.S. Federal law restricts those devices to sale by or on the order of a licensed dentist / physician.

We recommend an annual inspection of the prosthetic restoration by the dentist and the laboratory. This yearly inspection must consist of a screw check. If the screws are subject to unusual wear, the complete integrity of the implant abutment should be checked. New screws must be used for any adjustment correction, or replacement. Failure to follow this instruction puts the patient at risk and will void the warranty.

During any intraoral use and manipulation all DESS® products must be secured to prevent aspiration due to its shape and size.

INDICATIONS

Implant abutments are used for prosthetic restorations of dental implants or for assisting procedures in the dental laboratory.

CONTRAINDICATIONS

All materials used are biocompatible; however, some patients may be allergic or hypersensitive to any of the materials and their components (specified in the table).

COMPATIBILITY INFORMATION

All DESS® components are available for a variety of connections. For compatibility with dental implants and implant analogues please follow our catalogue or guides or contact your local distributor. To find necessary information about the ANKYLOS® C/X system, please read the last paragraph of these instructions.

STERILIZATION

All products are supplied non-sterile. For sterilization, autoclave treatment at 121 ° C for 30 minutes, drying 30 minutes (according to UNE-EN ISO 17665-1 and UNE-EN ISO 17665-2) is recommended.

RE-USE

Frequent use has a minimal effect on the instruments. The end of instrument life is usually determined by wear and tear and deterioration during use. Therefore, instruments can be reused with due care, provided they are not damaged or contaminated.

Discard single-use and worn reusable instruments immediately after use. Furthermore, any attempt to reuse such products seriously increases the risk of mechanical failure caused by material fatigue. No warranty claim resulting from the reuse of a single-use device will be accepted.

For the reuse of some instruments, they must be cleaned prior to sterilization. Here are some cleaning tips:

Never place instruments of different types of materials together.

Do not use metal brushes to remove impurities.

Use disposable syringes for instrument cavities.

When selecting detergents and disinfectants, make sure they are commercially available products and always follow the manufacturer's instructions.

SHELF LIFE

- ABUTMENTS/TIBASES:

The abutments and Ti-Base are intended to be handled and manipulated by professional dental technicians before being shaped into their final form. They are labelled "for single use", but do not degrade over time and are not susceptible to variations in environmental or storage conditions. Therefore, neither the shelf life of the device nor the shelf life of the package is limited and no "use by" date is provided on the label.

- SCREWS:

Clinical screws are intended to be handled by professionals for implant-prosthesis fixation, and once assembled, they are not susceptible to variations in environmental conditions. Their shelf life is not limited. Therefore, a date is not labelled.

| ANALOGUE / Class I | 14.XXX and 14.XXX-P10 | Surgical stainless steel AISI-303 Titanium alloy ELI Ti-6Al-4V | In the autoclave | Not recommended | ⊗ |
|--|---|--|-------------------------------------|--|---|
| TRANSFER / Class I | 17.XXX | Body: surgical stainless steel AISI-303 Screw: Titanium alloy ELI Ti-6Al-4V | In the autoclave | After sterilisation | - |
| SCREW / Class IIb | 19.XXX and 19.XXX-P10 | Titanium alloy ELI Ti-6Al-4V | Autoclave before use on the patient | Single use | ⊗ |
| CASTABLE ABUTMENT / Class IIb (28.XXX, 29.XXX), Class I (11.XXX-P10, 12.XXX-P10) | 28.XXX, 29.XXX and 11.XXX-P10, 12.XXX-P10, 33.XXX-P5, 32.XXX-P5 | Polyoxymethylene copolymer (POM) | N/A | Single use | ⊗ |
| Co-Cr BASE / Class IIb | 30.XXX, 31.XXX, 30.007H120, 30.007H122, 30.007H127, 30.007TB | Cobalt-Chromium-Molybdenum alloy | Autoclave before use on the patient | Single use | ⊗ |
| Ti BASE / Class IIb | 15.XXX, 16.XXX, 36.XXX, 37.XXX, 58.XXX, 59.XXX, 56.XXX, 57.XXX, 15.007H120, 15.007H122, 15.007H127, 15.007TB | Titanium alloy ELI Ti-6Al-4V | Autoclave before use on the patient | Single use | ⊗ |
| ABUTMENTS* (STRAIGHT AND ANGLED) / Class IIb | 13.XXX, 22.XXX and 23.XXX | Titanium alloy ELI Ti-6Al-4V | Autoclave before use on the patient | Single use | ⊗ |
| TRANSEPIHELIAL ABUTMENTS* (Multi-unit® and UniAbutment®) / Class IIb | 40.XXX/X, 41.XXX/X, 43.XXX/X and 43.XXX/XR | Titanium alloy ELI Ti-6Al-4V | Autoclave before use on the patient | Single use | ⊗ |
| HEALING ABUTMENT* / Class IIb | 20.XXX/X, 20.XXX, 20.7XX, 20.7XX/X, 20.100H120, 20.100H122, 20.100H127, 20.100TB | Titanium alloy ELI Ti-6Al-4V Polyether ether ketone (PEEK) | Autoclave before use on the patient | Single use | ⊗ |
| CONVERTER ABUTMENT / Class IIb | 34.001 and 35.001 | Titanium alloy ELI Ti-6Al-4V | Autoclave before use on the patient | Single use | ⊗ |
| TEMPORARY ABUTMENT* / Class IIb | 24.XXX, 25.XXX, 24.7XX, 25.7XX, 24.007H120, 24.007H122, 24.007H127, 24.007TB | Titanium alloy ELI Ti-6Al-4V Polyether ether ketone (PEEK) | Autoclave before use on the patient | Single use | ⊗ |
| CONOMETRIC ABUTMENT | 02.001_X 02.002_X | Titanium alloy ELI Ti-6Al-4V | Autoclave before use on the patient | Single use | ⊗ |
| SCANNING ABUTMENT / Class I | 50.XXX and 52.XXX, 53.XXX-P10 | Polyether ether ketone (PEEK) For ANK System: Ti alloy ELI Ti-6Al-4V Screw: Titanium alloy ELI Ti-6Al-4V | In the autoclave | In case of intraoral scanning: after sterilisation | - |
| DESSLOC** / Class IIb | 88.0XX/X and 88.XXX | Titanium alloy ELI Ti-6Al-4V, ZrN coated | Autoclave before use on the patient | Single use | ⊗ |
| DESSLOC® Retention Kit / Class I | DLRK-P1 and DLRK-P2, DLER-P1 and DLER-P2, DLN1-P4, DLN2-P4, DLN3-P4, DLN4-P4, DLN6-P4, DLN7-P4, DLRKL-P1 and DLRKL-P2, DLERL-P1 and DLERL-P2, DLN1L-P4, DLN2L-P4, DLN3L-P4, DLN4L-P4, DLN6L-P4, DLN7L-P4. | Insert - Polynil (Nylon) Cup - Titanium Alloy ELI Ti-6Al-4V Spacer ring - White Silicone | N/A | Single use | ⊗ |

| | | | | | |
|--------------------------------------|---|---|-------------------------------------|---|---|
| PRE-MILLED BLANKS / Class IIb | 61.XXX, 62.XXX, 63.XXX, 64.XXX 65.XXX, 66.XXX, and 71.XXX, 73.XXX | Titanium alloy ELI Ti-6Al-4V (61.XXX) and Cobalt-Chromium-Molybdenum alloy (71.XXX) | Autoclave before use on the patient | Single use | ⊗ |
| TOOLS / Class I | MDTXXX and DTXXXXXX | Handles: surgical stainless-steel SUS 316L Tips: stainless steel AISI 420 MOD | In the autoclave | In case of intraoral use: after sterilisation | - |
| TOOLS/Class IIa | DTXXXXXX, DTCXXXXXX and TTRXXXX | Torque application tools: AISI 420 MOD | In the autoclave | In case of intraoral use: after sterilisation | - |

APPLICATIONS AND DIRECTIONS FOR USE – ANALOGUES

The analogues are for simulating the dental implant connection and position on a stone master model. Before any laboratory manipulation ensure that the anti-rotation and anti-removal parts on the analogue are securely affixed. Check that the analogue and the prosthetic element connections match in size and type before tightening. A sealed and passive connection is needed. Do not re-use an analogue, because such a connection may be impossible to achieve with a product that may no longer conform to specifications.

APPLICATIONS AND DIRECTIONS FOR USE – TRANSFER

The transfer is used in the closed and open tray technique, for the transfer of the implant position from the intraoral situation to the model in the dental laboratory. Before being used, ensure that the connection seat for the implant is clean. Any traces of dirt could affect subsequent alignment of the prosthesis. Check compatibility with the implant model to which it will be connected. Remove the abutment, clean the connection with water/air and dry. Prepare the transfer in line with the implant connection and platform size. Place it on the implant and check that it is in the correct position. Put the screw in and tighten manually.

- For the closed tray technique choose the short screw. After tightening lock, the hex of the screw with wax. Use hydrocolloid, polyethylene, or soft silicone. After hardening place, the transfer in its groove and check the stability before sending to the laboratory. This technique is strongly recommended only for implants without axial divergence.
- For the open tray technique choose the long screw. After tightening position, the impression tray such as to ensure full access to transfer screws from the outside. Before taking impressions, the transfers can be ferulised with acrylic resin placed on dental floss between consecutive implants. After hardening release all screws and remove the impression tray.

APPLICATIONS AND DIRECTIONS FOR USE – CLINICAL SCREW

The screw is for fixing prosthetic restorations and auxiliary abutments over the implant.

or analogue. Make sure to secure the parts with their corresponding screws and observe the specified torque value placed on the label or ask your local DESS® distributor. You can download torque document from <https://www.dessdental.com/en-eu/downloads>. For best results, the following conditions must be adhered to meticulously:

- Use the suitable model key and size for tightening and unscrewing. If in doubt, check that the next size key does not fit into the seat. The screwdriver should be placed on the longitudinal axis of the prosthesis/implant assembly. A new screw should be used when assembling a prosthesis for the first time and for every check thereafter.
- For immediate load prostheses screw manually, avoiding excessive torque, and prevent the implant from turning while screwing.
- When transferring to the patient, do not use the same screw that was used in the laboratory.
- Make sure the correct model of screw is used for each case.
- Position the patient to avoid aspiration should the screw fall during screwing/unscrewing.**
- Check compatibility of the screw with the implant model to which it will be connected.

Some screw models are available with DLC coating [Diamond-Like Carbon, Chromium Carbonitride coating – CrCN] to provide a low friction surface, which improves the preload of the screw and provides better hold. The Chromium Carbonitride coating is non-toxic and has excellent mechanical, tribological and biological properties. DLC is often used in implants that make contact with blood such as stents and heart valves and to reduce wear in load bearing joints.

APPLICATIONS AND DIRECTIONS FOR USE – CASTABLE ABUTMEN

The castable abutment is for prosthetic restorations prepared by dental technicians in the dental laboratory. For greater connection precision use metal abutments which have been pre-milled. However, if you decide to use the casting method, please take the following precautions:

- Use wax that will create a sufficient gap around the abutment to compensate for the greater dilatation coefficient between the castable material and its surroundings.
- Use gentle torque to prevent deformation.
- Perform castings of a size and form that favours filling of the cavities and prevents the appearance of air bubbles.
- Use materials with high fluidity for castings in conflictive models.

APPLICATIONS AND DIRECTIONS FOR USE – Co-Cr BASE

The Co-Cr Base is used for prosthetic restorations prepared by dental technicians in the dental laboratory. The Co-Cr abutment is an implant abutment that consists of a cast-on Co-Cr alloy base and a plastic burnout sleeve. The laboratory processing of the Co-Cr base abutment uses the cast-on technique. During casting, the Co-Cr portion is suffused on the proposed contact areas by melting of the cast-on alloy. A metallic joint is obtained thanks to the retention of the grooves of the Co-Cr base and the suffused metal. **Processing:** The abutment must be carefully hand-tightened into the lab analogue with a lab screw. The lab screw is designed exclusively for laboratory processing of the abutment. The existing abutment screw should be used only for the final integration of the restoration after being sterilised. Create a wax collar of at least 0.3 mm in thickness above the base's metal indexing feature. It prevents cracking in the ceramic layer. Before placing casting sprues, the metal base should be cleaned with a cotton swab and alcohol. Casting sprues should be placed possibly in line with the longitudinal axis of the crown to avoid the bubble concentration inside the construction. The investment should freely pass through the screw channel. The elements need to be placed away from casting ring heating centre. While setting the sprue, ensure correct positioning of the wax- up casting object in the casting ring. The metal tube of the abutment should not be placed too close to the wall or floor of the ring to prevent heat removal (cooling vents). The use of speed investment is not recommended. The investment expansion should be adjusted to zero expansion. Make sure that the investment is free of air bubbles. Heat parameters need to be tuned for conventional CoCr casting to avoid any imperfections in the construction. Cooling the metal tube too rapidly can lead to defective casts. Observe the instructions of the investment material manufacturer. To ensure full casting, the final temperature during preheating of the casting muffle must be at least 950 °C. The final temperature should be held for 30 to 60 minutes according to the size of the casting ring. To prevent stresses in the cast structure resulting from cooling too rapidly, cool the casting muffle at room temperature. It is better to deflask the casting mould using an acid bath in an ultrasound cleaner. Sandblasting can damage the settlement. During neck polishing the connection surface needs to be covered with the analogue: this prevents any damage to the settlement part. If working with cemented restorations, you should develop crown and bridge frameworks on the abutments in the same way as for periodontic restorations.

COMPOSITION (percentage by mass)

Co 65.4 %
 Cr 27.75 %
 Mo 5.06 %
 Additional elements less than 1 % (Fe, Mn, Si, N, Ni, C, Ti, P, S)

TECHNICAL DATA

Yield strength¹ Rp 0.2 > 827 MPa
 Tensile strength¹ Rm > 1172 MPa
 Hardness¹ HV10 459
 Fracture elongation¹ A5 21 %
 Solidus¹ 1369 °C
 Liquidus¹ 1415 °C
 CTE¹ (25 – 500 °C) 14.1 x 10⁻⁶ K⁻¹
 Colour white

¹ Variations are possible depending on the batch.

APPLICATIONS AND DIRECTIONS FOR USE – TI-BASE

The Ti-Base is used for prosthetic restorations prepared by dental technicians in a dental laboratory. The main use of the titanium Ti-Base is to support the bridge or individual zirconium dioxide elements made by CAD/CAM systems or with manual milling machines. It can also be used as a connector between full-shape screwed crowns (from any material) and an implant, as a narrow abutment in some cases. **Contraindications:** do not use the titanium Ti-Base for restorations with cantilever on a single implant, bruxism, insufficient space or direct metal-to-Ti-Base casting.

To use this product in digital dentistry, the DESS® library and a scanning abutment are needed. To download and install one of our libraries specific to your dental CAD system please visit our website <https://www.dessdental.com/en-eu/downloads>

Some dental CAD/CAM systems require a double scan with a diagnostic wax-up. If you need additional assistance with installation and/or proper use of the files and with use of the piece itself, please contact your local DESS® distributor or find necessary information in corresponding instructions. All instructions can be downloaded from our website <https://www.dessdental.com/en-eu/downloads>. To improve cement adhesion, we recommend that it be thoroughly cleaned and degreased before cementation. You can use any implant luting cement (duo or chemical cured) extra-orally following the recommendations of the cement manufacturer. The internal surface of the zirconium structure (cementation zone) should be sandblasted and cleaned/degreased.

For secure grip the diameter and height of the Ti-Base should not be reduced.

APPLICATIONS AND DIRECTIONS FOR USE – STRAIGHT AND ANGLED CLINICAL ABUTMENTS

The straight and angled abutments are used for prosthetic restorations prepared by dental technicians in dental laboratory or by dentists chairside in dental clinics. Straight and angled abutments form the core of prosthetic crowns and bridges or are used for their support. Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts and that they are in the correct place and position for all secondary parts. We strongly recommend using radiography to verify the correct settlement after tightening. **Caution:** If customisation is required, do not grind the wall to less than 0.5 mm thick. **Warning:** Any grinding should be performed outside the mouth. If this is not possible, adequate suction must be used. Positioned the patient so that debris is not inhaled or swallowed. Check the torque according to implant connection and size. When cementing exclusively use luting cement or use temporary cement after sealing the access hole with Teflon tape or gutta-percha. Remove excess cement from edge of the crown to avoid peri-implantitis which can lead to implant loss. After cementation remember to check static and dynamic occlusion using the conventional protocol.

APPLICATIONS AND DIRECTIONS FOR USE – TRANSEPIHELIAL ABUTMENTS: Multi-unit® and UniAbutment®

These abutments are used for prosthetic restorations prepared by dentist in dental clinic. They are a pre-manufactured dental implant abutment directly connected to the dental implant intended for use as an aid in prosthetic rehabilitation for multiple fixed-removable or screw-retained prostheses. We strongly recommend checking the final abutment seating using radiographic imaging.

DESS® offers two systems for this kind of rehabilitation: Multi-unit® system and UniAbutment®. In the case of using the straight Multi-unit® abutment: place the appropriate abutment. Use plastic holder to facilitate the insertion. Tighten the abutment using DESS® key for Multi-unit® and manual torque wrench (for recommended torque use the Clinical Screw Recommended Torque document). In case of using the angled Multi-unit® abutment: place the appropriate abutment using metal holder at desired position and tighten the abutment using the system specific screwdriver in the manual torque wrench (for recommended torque use the Clinical Screw Recommended Torque document). Take an impression using the standard procedure and provisionalise the patient. If temporary prosthesis is not necessary, place the healing caps. In using the UniAbutment®: place the adequate abutment choosing the correct type (20° or 45° - depending on implant divergence) and height. Use the attached holder screw to facilitate insertion. Tighten the abutment using DESS® key for UniAbutment® and take impression using standard techniques (for recommended torque use the Clinical Screw Recommended Torque document).

APPLICATIONS AND DIRECTIONS FOR USE – HEALING ABUTMENT

The healing abutment is used on the patient as an auxiliary transmucosal abutment. It is placed on the implant fixture before prosthetic restoration to facilitate the formation of a soft tissue sulcus. Choose the correct height to ensure its correct function and prevent the transmission of chewing forces. Before placing ensure that implant platform stays free from any tissue residue. Use gentle manual torque when installing.

APPLICATIONS AND DIRECTIONS FOR USE – CONVERTER ABUTMENT

The converter abutment is used in prosthetic restorations prepared by dentists in dental clinics or by technicians in dental laboratories. Converter abutments are used for changing the connection from internal to external through the Morse taper connections to secure implant-to-abutment connection. These elements are placed over the dental implant acting as part of the prosthesis and are connected to the implant by a captive screw.

After placing in implant, use the correct DESS® part in a ratchet to tighten the abutment: for non-engaging abutments use an octagonal key. For engaging abutments use the corresponding screwdriver tip in adapter.

The recommended torque for both adapters can be found in the Clinical Screw Recommended Torque document, continue with standard prosthetic procedure.

APPLICATIONS AND DIRECTIONS FOR USE – TEMPORARY ABUTMENT

The temporary abutment is used for prosthetic restorations prepared by dentists at the chairside, it works as support for provisional screwed crowns, bridges and prosthesis. Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts. We strongly recommend using radiography to verify the correct settlement after tightening. **Warning:** Any grinding should be performed outside the mouth. If this is not possible, adequate suction must be used. Positioned the patient so that debris is not inhaled or swallowed. Use carbide burs and carborundum discs with water cooling. Before placing the temporary restoration check the torque according to implant connection and size. It is important to seal the access hole with Teflon tape or gutta-percha and composite. Any material excess near the implant platform should be removed, as it can provoke peri-implantitis and can lead to implant loss. It is very important to check static and dynamic occlusion avoiding excessive loading forces which may prolong or prevent the proper osseointegration, especially in cases of immediate load.

APPLICATIONS AND DIRECTIONS FOR USE – AURUMBase / ANGLEBase

The angled Ti-Base is used for prosthetic restorations prepared with axially tilted implants. The main use of the titanium Ti-Base is to correct this axis to support the bridge or individual zirconium dioxide elements fabricated in CAD/CAM or on the manual milling-copying machine. It can also be used as a connector between screw-retained full crowns (of any material) and an implant, or in some cases as a narrow clinical abutment. **Contraindications:** Do not use the titanium interface for single-unit restorations with extensions on a single implant, patients with bruxism, in cases of insufficient vertical space or for metal overcasts.

To use this product in digital dentistry you need the DESS® library and scanning abutments. To download and install one of our libraries specific to your CAD system please visit our website <https://www.dessdental.com/en-eu/downloads>

Some CAD/CAM systems require a double scan of the diagnostic wax-up. If you need further assistance with the installation and/or correct use of the libraries, or with the use of the part itself, please contact your local DESS® distributor. The necessary information can also be found in the corresponding instructions. All instructions can be downloaded from our website <https://www.dessdental.com/en-eu/downloads>

To improve the adhesion of the cement it is recommended to clean and degrease the surfaces prior to cementation. Any implant cement (dual-curing or self-adhesive) can be used for cementing on the model according to the cement manufacturer's instructions. The internal surface of the zirconium structure (cementation zone) area should be sandblasted and cleaned/degreased. For secure bonding, reduction of the interface diameter and height is not recommended.

APPLICATIONS AND DIRECTIONS FOR USE C-BASES

These attachments work just like straight Titanium Ti-Bases. They are used for prosthetic restorations performed in CAD/CAM with the CEREC and INLAB SW systems from Dentsply® Sirona. The function of these Ti-Bases is to support the cores of different materials (Zirconium, PMMA...) and to screw them to the implant. These restorations can also be anatomical crowns or anatomical bridges. To use this product in digital dentistry you have two options, use the DESS® library together with the DESS® scan-abutments or design with the Syrona® Cerec or Inlab software with your Syrona scan abutments. To download and install one of our libraries specific to your CAD system please visit our website <https://www.dessdental.com/en-eu/downloads>.

To improve the adhesion of the cement it is recommended to clean and degrease the surfaces before cementation. Any implant cement (dual-cure or self-adhesive) can be used for cementing onto the model according to the cement manufacturer's instructions. The ceramic surface in the cementation area should be sandblasted and cleaned/degreased. For safe bonding, reduction of the Ti-Bases diameter and height is not recommended.

APPLICATIONS AND DIRECTIONS FOR USE- ELLIPTIBase®

These are Ti-Bases designed for restorations with little interproximal space. This allows tighter restorations between the inferior-anterior group to be made without losing thickness of the nuclei or fully anatomical designs where they will be cemented in the ELLIPTIBase®. These Ti-Bases have two positions to improve their location on the implant. A (corner), position in the vertex of the hexagon B (Flat), position in the flat face of the hexagon. The function of these Ti-Bases is to support the nuclei of different materials (zirconium, PMMA...) and to screw them to the implants. These restorations can also be anatomical crowns or anatomical bridges. To use this product in digital dentistry it is not necessary to use the DESS® library. The CAD design of the elements to be manufactured will be carried out with the same software and library from Dentsply® Sirona.

To improve the adhesion of the cement it is recommended to clean and degrease the surfaces before cementation. Any implant cement (dual-cure or self-adhesive) can be used for cementing onto the model according to the cement manufacturer's instructions. The ceramic surface in the cementation area should be sandblasted and cleaned/degreased. For safe bonding, reduction of the Ti-

Bases diameter and height is not recommended.

APPLICATIONS AND DIRECTIONS FOR USE – SCANNING ABUTMENT

Scanning abutments are used for prosthetic restorations prepared by dental technicians in dental laboratories or by dentists at chairside in the CAD/CAM process. Designed to obtain geometric data from the master model using a desktop laboratory 3D scanner or for optical impressions using an intraoral 3D scanner. For proper processing a DESS® digital library is needed. To download and install one of our libraries specific to your dental CAD system please visit our website <https://www.dessdental.com/en-eu/downloads>. Before scanning the abutment and analogue they should be visually checked for surface damage or implant connection for any tissue residue. For greater scanning precision, we recommend locating the flat surface of the scanning abutment in palatal/lingual orientation. Fasten the abutment using the corresponding screw (directions in our catalogue and on the product label) by hand or with a maximum torque of 10 Ncm. In case of intraoral scan abutment, fix it with fingers using the embedded screw. The scanning abutment is a precision tool and being overtightened may change its geometry causing errors in the scanning process and discrepancy in accuracy. For most of scanners a scan-enhancing spray is not required. The scanning process should be performed as recommended by the CAD/CAM system manufacturer. It is important to choose the correct implant connection in the software and the corresponding type for the chosen restoration (engaging/non-engaging). After the process the scanning abutment can be loosened and placed gently on the tray or box. To complete the scan, gathering additional information (e.g., occlusion index, silicone bite, gingiva shape) may be required. In case of intraoral use sterilisation is required.

APPLICATIONS AND DIRECTIONS FOR USE – DESSLoc®

DESSLoc® is a resilient implant attachment used for prosthetic restorations prepared by dental technicians in dental laboratories or by dentist at chairside. Recommended for removable dentures. **Procedure:** Measure the tissue thickness from the implant platform to the highest side of the gingival crest near the implant. Choose the corresponding height or 1.5 mm higher (the attachment should not be submerged) and put in place using the recommended torque from the Clinical Screw Recommended Torque document. The retention male part can be mounted in the dental clinic as well as in dental laboratory with self-curing resin using standard techniques.

APPLICATIONS AND DIRECTIONS FOR USE – DESSLoc® RETENTION KIT

DESSLoc® retention system is a part of the anchoring system for full or partial overdentures on implants, in cases where conventional full dentures have problems with retention and stability. The male retention kit can be mounted in the dental clinic using the self-curing resin, and in the dental laboratory, using the standard procedure. **Direct method:** Place the silicone spacer over the attachment covering the entire area around it. If there is any space between the cup and soft tissue, put more spacers until the space disappears. Place the cup with the black Nylon male part on the DESSLoc® abutment in the patient's mouth. Prepare a site for the cup in the prosthesis and perforate the buccal side to evacuate excess material. Apply the dual- or cold-curing resin around the cup and within the holes. Place the prosthesis in place while keeping the height of the occlusion and hold in position until the resin hardens. Remove the denture, spacers and cut off the excess. Polish before changing to the final male Nylon. **Indirect method:** same as direct, but on the master model with 14.054 analogues. **Retention force/divergence between implants:** blue (1.5 lbs/680g, up to 20°), pink (3 lbs/1360g, up to 20°), transparent (5 lbs/2270g, up to 20°), black (Processing Nylon male, 0 lbs, 0g, up to 40°), red (1 lbs/450g, up to 40°), orange (2 lbs/910g, up to 40°), green (4 lbs/1810g, up to 40°). DLRK-P1 and DLRK-P2 are kits up to 20°, DLER-P1 and DLER-P2 are kits up to 40°.

APPLICATIONS AND DIRECTIONS FOR USE – MASSIVE ABUTMENT

Massive abutments are used for prosthetic restorations prepared by dental technicians in dental laboratories. The pre-milled implant-abutment Ti-Base of the blank is fabricated on an industrial lathe in order to provide consistent and accurate results. The connection is more accurate, and the abutment shape can be milled with burs to meet specific patient requirements.

Caution: For proper milling a DESS® fixture adapter is needed. The minimum wall thickness shall not be less than 0.45 mm.

APPLICATIONS AND DIRECTIONS FOR USE – TOOLS

Torque transmission instruments are used to apply or transmit torque to instruments, implants or prosthetic devices, which can be used in conjunction with a ratchet or surgical handpiece. All instruments used intraorally should be sterilized by the practitioner prior to use.


Torque transmission instruments such as interchangeable tips can be used mechanically or manually.

If used manually, it must be ensured that the interchangeable tips are securely locked in the handle. The handle has a groove for floss clamping to prevent aspiration by the patient.

When using the interchangeable tips, make sure that they match in size and shape to the shape of the screw head. Periodic replacement of this type of instruments is recommended. Before using any instruments, check that they are in perfect condition and suitable for use.

In the case of adapters, check if it corresponds to the 4x4 connection and its system.

(multi-unite®, UniAbutment®, Octagonal or Locator®).

INTERNAL ANK SPECIAL REQUIREMENTS

The thickness of the anti-rotation elements is reduced due to the widened screw channel in order to accommodate the use of a normal screw. To prevent bending of those parts under applied pressure, please ensure that the screw is threaded through the abutment whilst restoration takes place and before placing it in the analogue or implant. If this precaution is taken, the anti-rotation elements remain supported and safe from shear and compression forces. If, for any reason the screw needs to be removed, place it back again in its position before applying any loads on the abutment.

DESS® MILLING TOOL

The milling tool can create bars, bridges and multiples independently. For those professionals who wish to create their own restorations, the DESS® milling tool allows you to create angled channel solutions up to 25°. The professional can carry out their own project and ensure the required precision. This new DESS® tool accurately reproduces the platform of the Torx screw for AURUMBase®, the DESS® angled channel solution. Below are the compatibilities VHF®:3MM, STRAUMANN®: 3MM, IMES-ICORE®:3MM & 6MM, ZIRKONZAHN®:3MM 6MM, ROLAND®: 4MM, YENADENT®: 4MM & 6 MM.

NOTICE ON SERIOUS INCIDENTS

For patients / users / third parties in the European Union with an identical regulatory regime (Regulation 2017/745/EU) if, due to the use of the product, a serious incident occurs, notify the manufacturer and your national authority. The manufacturer's contact information for the notification of this type of incident is:

Terrats Medical SL
Calle Mogoda 75-99
08210 – Barberá del Vallés -Barcelona (Spain)
Tel + 34 935 646 006
cesar.escribano@dessdental.com

EXPLANATION OF THE SYMBOLS THAT APPEAR ON LABELS

| | | |
|---|--|----------------------|
|  For single use only | LOT Batch No. | REF Reference |
|  NON Sterile product |  Read instructions | CE 0051 CE |
| Rx Only Use only under prescription | Medical product compliance marking controlled by the notified body number 0051 (Class IIb and Class IIa) | |
|  Manufacturer information | Manufacturer: Terrats Medical SL Calle Mogoda 75-99 08210 – Barberá del Vallés -Barcelona (Spain) Tel + 34 935 646 006. Fax +34 935 647 317 info@dessdental.com | |
| CH REP | Innovdentaire Route des jeunes 41A, 1227 CAROUGE, Switzerland regis@innovdentaire.ch | |

All DESS® products are manufactured in accordance with ISO 9001 and ISO 13485 and bare the CE marking. 3i® and 3i® Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations, USA
 ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XIVE® and FRIADENT® and UniAbutment® are registered trademarks of DENTSPLY Implants
 BioHorizons® is a registered trademark of BioHorizons, INC.
 Brånemark System®, NobelReplace® and Replace Select™ and NobelActive™ and NobelReplace® CC and NobelSpeedy™ and Multi-unit® are registered trademarks of Nobel Biocare AB, Sweden
 ITI®, Straumann® and synOcta® are registered trademarks of Straumann Holding AG, Switzerland
 MIS® is a registered trademark of MIS Implant Technologies, Ltd
 Zimmer Screw-Vent® and SwissPlus® are registered trademarks of Zimmer Dental INC, USA
 Dyna® is a registered trademark of Dyna Dental Engineering BV
 Locator® is a registered trademark of Zest Anchors INC.

The photographs of the products are for information purposes only. The type, reference and connection of each element should be checked. Some DESS® items may not be available in your country. These instructions replace all previous editions. For any detailed information about DESS® products please contact your local distributor.

Medical product compliance marking controlled by the notified body number 0051 (Class IIb and class IIa)

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